

Accessories

For Atrium Thoracic Chest Drainage Systems



19910 – Red cap

- For capping patient line before disposal
- For use if only one patient line is needed on a dual drain
- Twenty-five (25) each per case – sterile



19913 – Adult / pediatric connector 1/4" x 3/8"

- For use with an adult drain and smaller catheters – 16, 20 French
- Twenty-five (25) each per case – sterile



19914 – Standard adult connector 3/8" x 3/8"

- Standard 5-in-1 connector
- For use on larger catheters – 24 French and larger
- Twenty-five (25) each per case – sterile




19915 – Pediatric connector 1/4" x 1/4"

- Larger connector included with infant / pediatric drain
- For use with an adult drain and smaller catheters – 16 French or larger thoracic catheters, 9 and 11 mm mediastinal catheters
- Twenty-five (25) each per case – sterile



19916 – Sims connector 3/16" x 1/4"

- Same connector included with infant / pediatric drain
- For use with smaller catheter – 8, 10, 12 French
- Twenty-five (25) each per case – sterile



19917 – Y connector 3/8" x 3/8" x 3/8"

- Connects two chest tubes to one adult chest drain
- Twenty-five (25) each per case – sterile
- For use with 11 mm mediastinal catheter

Accessories



19918 – Step down connector 3/16" x 3/8"

- For use with adult chest drain
- For use when connecting Blake drain to adult drain
- Fits smaller catheters – 8, 10, 12 French



19907 – Patient tube pack

- 3/8" ID tubing
- Safety stock if patient tube sterility is accidentally compromised
- Provides cost savings to patient and hospital
- Not made with natural rubber latex
- Ten (10) each per case – sterile



19909 – Y extension set

- Simple, pre-assembled extension connection for attaching two chest tubes to one chest drain
- 3/8" ID tubing
- Connectors fit 24 French and larger
- Ten (10) each per case – sterile



19927 – Luer-lock connector with tube

- Allows connection of a standard chest drain or mobile drain to a luer-lock catheter
- Twenty-five (25) each per case – sterile



19928 – Suction Y adapter

- For use when connecting two chest drains to one suction source
- Ten (10) each per case – sterile



19930 – Express suction bulb

- For use with the Express Mini 500 (No. 16400)
- Allows portable, manual vacuum application
- Ten (10) each per case – sterile

This document is intended to provide information to an international audience outside of the US.

Manufacturer · Atrium Medical Corporation · 40 Continental Blvd · Merrimack, NH 03054 · USA

www.getinge.com

GETINGE 

Protected by the following international and U.S. patent(s): patent(s): patents.maquet.com · CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician. · Getinge Group, Getinge, Arjo-Huntleigh, Atrium and Maquet are trademarks or registered trademarks of Getinge AB, its subsidiaries or affiliates in the United States or other countries · Getinge, Arjo, Atrium and Maquet are registered with the U.S. Patent and Trademark Office · Copyright 2014 Atrium Medical Corp. or its affiliates · All rights not expressly granted are reserved · Refer to instructions for use for current indications, warnings, contraindications, and precautions · Part No. 010112 Rev AC · Refer to Instructions for Use for current indications, warnings, contraindications, and precautions.